## Histoplasma Urine Antigen LFA Control Kit

### 1. Intended Use

The Histoplasma Urine Antigen LFA Assay Controls are quality control reagents for use only with the Histoplasma Urine Antigen LFA Test Kit.

### 2. Reagents

Histoplasma Urine Antigen LFA Control Kit, 3 tests, Catalog # 1020

Store the kit at 2°C to 8°C.

Component	Contents	Quantity
Negative Control	Non-human derived; Synthetic urine	1 X 0.330 mL
Positive Control	Non-human derived galactomannan in synthetic urine matrix	1 X 0.330 mL

### 3. Materials Required But Not Provided

- A. Histoplasma Urine Antigen LFA Test Kit.
- B. Proper Personal Protective Equipment (PPE), including lab coat, gloves, and eye protection.

### 4. Warnings and Precautions

- A. For *in vitro* diagnostic use only. For professional use only.
- B. The package inserts for both the Histoplasma Urine Antigen LFA Test Kit and the Histoplasma Urine Antigen LFA Control Kit should be read completely prior to testing. The user must adhere to the method instructions or risk adversely affecting the assay performance characteristics.
- C. When handling the Controls, wear a lab coat, protective gloves, and eye protection. Use universal standard precautions when handling the kit. Wash hands thoroughly after using the kit.
- D. Do not eat, drink, or smoke in any area where the kit is being used.
- E. Biological spills should be wiped with a disinfectant, such as 10% bleach or 70% ethanol. Materials used to wipe the spills may require biological waste disposal.
- F. Dispose of all materials used in running the assay as if they are biohazardous. Handle and dispose of laboratory chemical and biological waste in accordance with applicable laws and regulations.
- G. Do not use any components in the kit that are beyond their expiry.
- H. Cross contamination can occur in a negative assay due to improper handling or incorrect pipetting.

## 5. Reagent Handling and Storage

- A. Store Controls at 2°C to 8°C.
- B. Allow the Controls to come to room temperature before use.
- C. Controls are acceptable for use until the expiry listed on the package label.
- D. Sufficient samples are provided to test the Controls 3 times.

## 6. Procedure

Refer to the Histoplasma Urine Antigen LFA Test Kit package insert for procedural instructions.

## 7. Assay Interpretation

- A. Assay Acceptance:
  - 1. Determine whether there is a visible test line at the “T” location in the result window and a visible control line at the “C” location in the result window.
  - 2. Absence of a control line indicates an INVALID result. The assay should be repeated.
  - 3. If the sample fails to flow properly or if the assay exhibits anomalies such as excessively high background, the test is invalid and the test line is not evaluated. The assay should be repeated.
- B. Expected Results:
  - 1. Negative Control: Absence of a test line and presence of a control line.
  - 2. Positive Control: Presence of a test line and presence of a control line. If the test line is present but appears uneven or spotty, it may still be indicative of a positive test line.

## 8. Limitations of the Procedure

- A. The Histoplasma Urine Antigen LFA Control Kit is for use only with the Histoplasma Urine Antigen LFA Test Kit.
- B. Color blindness may interfere with the ability to interpret results.

	Manufacturer	<b>IVD</b>	In Vitro Diagnostic Device
	Expiration Date	<b>REF</b>	Reference Number
	Storage Temperature	<b>LOT</b>	Lot Number
	Consult Instructions for Use	<b>CE</b>	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices
	Single Use. Do Not Reuse		Do Not Use if Package is Damaged
<b>EC</b> <b>REP</b>	Authorized Representative in the European Community		

MiraVista Diagnostics, LLC  
 4705 Decatur Blvd.  
 Indianapolis IN 46241 USA  
 (866) 647-2847  
[www.miravistalabs.com](http://www.miravistalabs.com)

CEpartner4U BV  
 Esdoornlaan 13  
 3951 DB Maarn  
 The Netherlands  
[www.cepartner4u.com](http://www.cepartner4u.com)

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