Histoplasma Urine Antigen LFA Test Kit

1. Intended Use

The Histoplasma Urine Antigen LFA is an immunochromatographic lateral flow assay for the qualitative detection of *Histoplasma* antigen in urine samples.

2. Indications for Use

The Histoplasma Urine Antigen LFA can be used to provide presumptive evidence of infection with *Histoplasma* and aid in the diagnosis of histoplasmosis in persons with signs and symptoms of histoplasmosis.

3. Summary and Explanation

Histoplasmosis is an infection caused by the dimorphic fungus *Histoplasma capsulatum*. *Histoplasma capsulatum* exists throughout the world but is highly endemic to certain parts of North America, Central, and South America, and much less so in parts of Asia, Southeast Asia, China, India, Australia, Africa and Europe. In the United States, endemic regions include the Ohio River Valley and lower Mississippi River Valley¹. Infection mainly occurs through inhalation of aerosolized microconidia (spores). Infection usually starts in the pulmonary region but may disseminate to other regions of the body in more serious cases. Immunocompromised individuals, such as those with AIDS or those who have recently undergone organ or bone marrow transplant, are more susceptible to disseminated infection. Disseminated histoplasmosis, if left untreated, may be fatal². Diagnosis of histoplasmosis is usually dependent on multiple factors, including clinical symptoms, serology, histopathology, antigen detection, culture, and chest x-rays^{1,3,4}.

4. Principle of the Procedure

The assay is a lateral flow-based immunoassay which employs polyclonal antibodies in the direct detection of *Histoplasma* antigen in urine. Urine sample is added to Sample Diluent using the empty tube provided in the kit. The diluted sample is then added to the sample port of the device. If *Histoplasma* antigen is present in the diluted sample, the antigen binds to an anti-*Histoplasma* polyclonal antibody-colloidal gold conjugate as the sample moves through the device. The *Histoplasma*-capture antibody bound to the assay membrane at the Test position of the device reading window binds the antigen-antibody-colloidal gold complex and yields a visible line. As an assay control, anti-rabbit antibody bound at the Control position of the reading window yields a visible line, indicating adequate flow of diluted sample through the test device.

5. Reagents

Histoplasma Urine Antigen LFA Test Kit, 25 tests, Catalog # 1010
 Store the kit at room temperature (15°C to 30°C).

Component	Contents	Quantity
<i>Histoplasma</i> Antigen Lateral Flow Device	An immunochromatographic assay device for the detection of <i>Histoplasma</i> antigen in urine Sealed in a pouch with desiccant pillow	1 device per pouch, 25 devices per kit
Sample Diluent	Urine Sample Diluent	1 X 1.0 mL
Tube	1.5 mL Tube	25
Package Insert	Instructions for Use and Warnings and Precautions	1

6. Materials Required But Not Provided

- A. Proper Personal Protective Equipment (PPE), including lab coat, gloves, and eye protection.
- B. Calibrated single-channel, variable volume micropipette(s).

7. Warnings and Precautions

- A. For *in vitro* diagnostic use only. For professional use only.
- B. The package insert should be read completely prior to testing samples.
- C. All urine samples should be handled as if they are capable of transmitting infection.
- D. Urine is the only matrix acceptable for testing using this assay.
- E. When handling any of the kit components, wear a lab coat, protective gloves, and eye protection.
- F. Patient samples may contain infectious agents. Patient specimens and all materials which they contact should be handled at Biosafety Level 2 as recommended in the CDC/NIH manual “Biosafety in Microbiology and Biomedical Laboratories.”
- G. Wash hands thoroughly after using the kit.
- H. Do not eat, drink, or smoke in any area where the kit is being used or where samples are being handled.
- I. Biological spills should be disinfected using solutions such as 10% bleach or 70% ethanol. Materials used to wipe the spills may require biological waste disposal.
- J. Dispose of all urine samples and materials used in running the assay as if they are biohazardous. Handle and dispose of laboratory chemical and biological waste in accordance with applicable laws and regulations.
- K. Refer to the safety data sheet (SDS) for more information.
- L. Inspect lateral flow device pouch for any defects prior to use. Do not use device if the pouch has any tears or holes.
- M. Do not use any components in the kit that are beyond their expiry.
- N. Use care when pipetting samples. Failure to properly add samples to Sample Diluent or to the device may increase the chance of false negatives or false positives. If a procedural error is suspected, repeat testing should be performed.
- O. Cross contamination can occur due to improper handling or incorrect pipetting.

8. Reagent Handling and Storage

- A. Store all kit reagents at room temperature (15°C to 30°C).
- B. LFA kits are acceptable for use until the expiry listed on the kit label.
- C. Use lateral flow device within 1 hour after removal from pouch.

9. Urine Sample Handling and Storage

- A. Urine is the only sample type acceptable for use with the *Histoplasma* Urine Antigen LFA Test Kit.
- B. Store urine samples at 2°C to 8°C for up to 14 days or frozen (-20°C to -80°C) if the sample is to be stored longer. **Urine samples that are improperly stored may give false positive results as they may contain overgrowth of fungus and/or bacteria.**
- C. Avoid the use of sample storage containers that contain the preservatives boric acid, sodium borate, or sodium formate, such as culture and sensitivity preservative tubes.
- D. Test urine samples within 2 hours after the addition of Sample Diluent.

10. Procedure

- A. Use this procedure for the testing of urine samples and Assay Controls.
- B. Bring urine sample(s) to room temperature (15°C to 30°C).
- C. Remove the lateral flow device from the foil pouch.
- D. Using a micropipette, transfer 25 µL of Sample Diluent to a provided tube.
- E. Using a micropipette with a new pipette tip, transfer 100 µL of urine sample to the tube of Sample Diluent.
- F. Mix diluted sample tube gently for 5 seconds. Ensure urine sample and Sample Diluent are completely mixed.
- G. Using a micropipette with a new pipette tip, add 100 µL of diluted urine sample to the lateral flow device sample reservoir.
- H. Allow assay to process at room temperature (15°C to 30°C) for 40 ± 10 minutes.

11. Assay Interpretation

- A. Assay Acceptance:
 - 1. Determine whether there is a visible pink test line at the “T” location in the result window and a visible red control line at the “C” location in the result window.
 - 2. Absence of a control line indicates an INVALID result. The assay should be repeated.
 - 3. If the sample fails to flow properly or if the assay exhibits anomalies such as excessively high background, the test is invalid and the test line is not evaluated. The assay should be repeated.
- B. Sample Interpretation:
 - 1. Absence of a test line and presence of a control line is interpreted as NEGATIVE.
 - 2. Presence of a test line and presence of a control line is interpreted as POSITIVE. If the test line is present but appears uneven or spotty, it may still be indicative of a positive test line.

12. Limitations of the Procedure

- A. A negative result does not rule out the possibility of the diagnosis of histoplasmosis.
- B. The Histoplasma Urine Antigen LFA Test Kit has not been evaluated for use with sample types other than urine.
- C. The antibody in this device has been shown to cross-react with *Aspergillus* spp., *Blastomyces dermatitidis*, *Coccidioides* spp., *Cryptococcus neoformans*, Cytomegalovirus, *Enterococcus faecalis*, *Histoplasma duboisii*, *Klebsiella pneumoniae*, *Paracoccidioides brasiliensis*, *Sporothrix schenckii*, *Staphylococcus aureus*, *Streptococcus agalactiae*, and *Talaromyces marneffeii* and may result in a positive test.
- D. It is possible that increased levels of urine protein can interfere with the Histoplasma Urine Antigen LFA and cause a false positive result.
- E. Above normal levels of the following substances may interfere with Histoplasma Urine Antigen LFA sensitivity: bilirubin, creatinine, glucose, ketones, and urobilinogen.
- F. Color blindness may interfere with the ability to interpret results.

13. Performance Data

A. Precision and Repeatability

i. Precision

Within-lab precision was determined by testing 4 test urine samples—a moderate positive, a low positive, a high negative, and a negative—in duplicate over 10 runs, with two runs per day for 5 days. 100% of samples were correctly interpreted. See Table 13.1.

Table 13.1. Within-lab Precision Results.

Sample	Expected Result	Day 1		Day 2		Day 3		Day 4		Day 5	
		AM Run	PM Run	AM Run	PM Run	AM Run	PM Run	AM Run	PM Run	AM Run	PM Run
Negative Control	Negative	NEG	NEG	NEG	NEG	NEG	NEG	NEG	NEG	NEG	NEG
Positive Control	Positive	POS	POS	POS	POS	POS	POS	POS	POS	POS	POS
Mod Pos 1	Positive	POS	POS	POS	POS	POS	POS	POS	POS	POS	POS
Mod Pos 2	Positive	POS	POS	POS	POS	POS	POS	POS	POS	POS	POS
Low Pos 1	Positive	POS	POS	POS	POS	POS	POS	POS	POS	POS	POS
Low Pos 2	Positive	POS	POS	POS	POS	POS	POS	POS	POS	POS	POS
High Neg 1	Negative	NEG	NEG	NEG	NEG	NEG	NEG	NEG	NEG	NEG	NEG
High Neg 2	Negative	NEG	NEG	NEG	NEG	NEG	NEG	NEG	NEG	NEG	NEG
Neg 1	Negative	NEG	NEG	NEG	NEG	NEG	NEG	NEG	NEG	NEG	NEG
Neg 2	Negative	NEG	NEG	NEG	NEG	NEG	NEG	NEG	NEG	NEG	NEG
Percent Correlation		100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0

ii. Repeatability

Within-run repeatability was determined by testing 20 replicates of a low positive sample and 20 replicates of a negative urine sample in a single run. 100% of replicates were correctly interpreted.

B. Analytical Sensitivity

i. Limit of Detection

In order to determine LFA analytical sensitivity, urine samples containing 0.6 ng/ mL – 2.2 ng/ mL of *Histoplasma* antigen were tested in 40 replicates each in each of 3 lots of Histoplasma Antigen lateral flow devices. The lowest concentration at which at least 95% of the replicates were interpreted as positive was selected as the assay limit of detection (LoD) for each lot. The highest of the three observed lot limits of detection was selected as the assay LoD. The LoD was determined to be 1.8 ng/ mL.

ii. Interfering Substances

The presence of substances that could potentially interfere with assay sensitivity was investigated by testing abnormal urines spiked with *Histoplasma* antigen. See Table 13.2.

Table 13.2. Interfering substances sensitivity study.

Substance	# Negative / Total
Antihistamine	0/12
Ascorbic Acid	0/12
Bilirubin	3/12
Blood/ Hemoglobin	0/12
BUN	0/12
Cough Syrup	0/12
Creatinine	4/12
Glucose	10/14
Ketones	2/10
Microalbumin	0/12
Nitrite	0/10
Protein	0/16
Uric Acid	0/12
Urobilinogen	3/10
WBC	0/10

C. Analytical Specificity

i. Interfering Substances

The presence of substances that could potentially interfere with assay specificity was investigated by testing abnormal urines. See Table 13.3.

Table 13.3. Interfering substances specificity study.

Substance	# Positive / Total
Antihistamine	0/4
Ascorbic Acid	0/4
Bilirubin	0/6
Blood/ Hemoglobin	0/4
BUN	0/4
Cough Syrup	0/4
Creatinine	0/6
Glucose	0/7
Ketones	0/5
Microalbumin	0/6
Nitrite	0/5
Protein	1/8
Uric Acid	0/6
Urobilinogen	0/4
WBC	0/4

ii. Cross-Reacting Organisms

The effect of potentially cross-reacting organisms on assay specificity was investigated. See Table 13.4.

Table 13.4. Cross-Reacting Organisms study.

Organism	# Positive/ Total
<i>Aspergillus flavus</i>	2/2
<i>Aspergillus fumigatis</i>	2/2
<i>Aspergillus nidulans</i>	1/2
<i>Aspergillus niger</i>	2/2
<i>Aspergillus terreus</i>	2/2
<i>Blastomyces dermatitidis</i>	2/2
<i>Candida albicans</i>	0/2
<i>Candida glabrata</i>	0/2
<i>Candida krusei</i>	0/2
<i>Candida parapsilosis</i>	0/2
<i>Candida tropicalis</i>	0/2
<i>Coccidioides immitis and posadasii</i>	2/2

Organism	# Positive/ Total
<i>Cryptococcus neoformans</i>	1/2
Cytomegalovirus	1/2
<i>Enterococcus faecalis</i>	1/2
<i>Enterococcus faecium</i>	0/2
<i>Escherichia coli</i>	0/2
<i>Histoplasma duboisii</i>	2/2
<i>Klebsiella pneumoniae</i>	1/2
<i>Paracoccidioides brasiliensis</i>	2/2
<i>Proteus mirabilis</i>	0/2
<i>Pseudomonas aeruginosa</i>	0/2
<i>Sporothrix schenckii</i>	2/2
<i>Staphylococcus aureus</i>	2/2
<i>Staphylococcus saprophyticus</i>	0/2
<i>Streptococcus agalactiae</i>	1/2
<i>Talaromyces marneffeii</i>	2/2

D. Clinical Sensitivity and Specificity

To establish clinical sensitivity and specificity, urine samples were tested from three endemic sites in the United States. Patients with a culture positive for *Histoplasma* served as cases (n = 44), while patients with a culture negative for *Histoplasma* served as controls (n = 296). One urine sample was tested from each patient. Histoplasma Urine Antigen LFA results were compared to culture as the reference standard. See Tables 13.5 and 13.6.

Table 13.5. Histoplasma Urine Antigen LFA clinical sensitivity and specificity diagnostic table.

<i>Histoplasma</i> Urine Antigen LFA	Culture	
	Positive (Case)	Negative (Control)
Positive	41	9
Negative	3	287
TOTAL	44	296

Table 13.6. Histoplasma Urine Antigen LFA clinical sensitivity and specificity with 95% confidence intervals.

Sensitivity	93.18% (81.34% to 98.57%)
Specificity	96.96% (94.31% to 98.60%)
Positive Likelihood Ratio	30.65 (16.03 to 58.60)
Negative Likelihood Ratio	0.07 (0.02 to 0.21)
Disease Prevalence	12.94% (9.56% to 16.98%)
Positive Predictive Value	82.00% (68.56% to 91.42%)
Negative Predictive Value	98.97% (97.01% to 99.79%)

14. References

- (1) Kauffman CA. Histoplasmosis: a clinical and laboratory update. *Clin Microbiol Rev* **2007 Jan**;20(1):115-32.
- (2) Wheat LJ, Slama TG, Norton JA, et al. Risk factors for disseminated or fatal histoplasmosis. Analysis of a large urban outbreak. *Ann Intern Med* **1982 Feb**;96(2):159-63.
- (3) Hoz RL, Loyd JE, Wheat LJ, Baddley JW. How I Treat Histoplasmosis. *Curr Fungal Infect Rep* **2013 Mar**;7(1):36-43.
- (4) Hage CA, Ribes JA, Wengenack NL, et al. A multicenter evaluation of tests for diagnosis of histoplasmosis. *Clin Infect Dis* **2011 Sep**;53(5):448-54.

	Manufacturer	IVD	In Vitro Diagnostic Device
	Expiration Date	REF	Reference Number
	Storage Temperature	LOT	Lot Number
	Consult Instructions for Use	CE	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices
	Single Use. Do Not Reuse		Do Not Use if Package is Damaged
EC REP	Authorized Representative in the European Community		

MiraVista Diagnostics, LLC
 4705 Decatur Blvd.
 Indianapolis IN 46241 USA
 (866) 647-2847
www.miravistalabs.com

CEpartner4U BV
 Esdoornlaan 13
 3951 DB Maarn
 The Netherlands
www.cepartner4u.com

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14. Referencias

- (1) Kauffman CA. Histoplasmosis: a clinical and laboratory update. Clin Microbiol Rev **2007** Jan;20(1):115-32.
- (2) Wheat LJ, Slama TG, Norton JA, et al. Risk factors for disseminated or fatal histoplasmosis. Analysis of a large urban outbreak. Ann Intern Med **1982** Feb;96(2):159-63.
- (3) Hoz RL, Loyd JE, Wheat LJ, Baddley JW. How I Treat Histoplasmosis. Curr Fungal Infect Rep **2013** Mar;7(1):36-43.
- (4) Hage CA, Ribes JA, Wengenack NL, et al. A multicenter evaluation of tests for diagnosis of histoplasmosis. Clin Infect Dis **2011** Sep;53(5):448-54.

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